

TRAUMA OXYGEN UNIT

SPECIFICATIONS

- **Unit Weight:** Approximately 20 pounds
- **Unit Size:** 8 1/2" x 8 1/2" x 23"
- **Cylinder Capacity:** Approximately 19 cu. ft. of U.S.P. oxygen (529 liters) at 2015 psi
- **Unit Case:** Cylinder and Regulator are protected by a shock-resistant portable polyethylene case. Provides space for storage of first aid supplies and personal/clean-up equipment.
- **Can be wall mounted or cart wheeled for ease of mobility.**



For more information contact LifeSafe Services at:

LifeSafe Services, Inc.
5971 Powers Ave., Suite 8
Jacksonville, FL 32217 USA
Phone: 1.904.730.4800 • 1.888.767.0050
Fax: 1.904.730.4801 • 1.888.767.4018
Email: terry.timmins@lifesafeservices.com
www.lifesafeservices.com