

## ERIE PORTABLE OXYGEN UNIT

### SPECIFICATIONS

- Filled and ready to use with 246 liters of medical oxygen USP – over a 40 minute supply.
- Permanently attached on/off handle – no wrench needed.
- 6 liter fix flow regulator is standard. Can also be fitted with a variable flow regulator.
- Pressure gauge in cylinder post valve – shows oxygen level in cylinder.
- Medium concentration mask and 7' kink proof hose – Sanitary single patient use.
- Refillable light weight aluminum cylinder – Maximum economy.
- Highly visible instructions.
- Complete system (246 liter unit) weighs only 7 pounds.
- Meets FDA Requirements – Use without a prescription.
- Durable shoulder bag which leaves hands free or durable case for easy wall mounting available at additional charge.



For more information contact LifeSafe Services at:

**LifeSafe Services, Inc.**  
5971 Powers Ave., Suite 8  
Jacksonville, FL 32217 USA  
Phone: 1.904.730.4800 • 1.888.767.0050  
Fax: 1.904.730.4801 • 1.888.767.4018  
Email: [terry.timmins@lifesafeservices.com](mailto:terry.timmins@lifesafeservices.com)  
[www.lifesafeservices.com](http://www.lifesafeservices.com)